



# Obstetrics-Pelvimeters



**04-6379**  
**THOLE**  
for internal and external measurements, graduated in centimeters

1/2



**04-6380**  
**DE LEE**  
graduated in centimeters and inches

1/2



**04-6381**  
**COLLIN**  
graduated in centimeters

1/2



**04-6382**  
**MARTIN**  
graduated in centimeters

1/2



1/2



**04-6383**  
**SCHNEIDER**  
graduated in centimeters

for internal measurements between symphysis and promontory

The additional straight arm can be replaced. The pelvimeter with straight arm is used for localisation of the portio vaginalis uteri and of vaginal metastasis.




**04-6384**  
**BREISKY**  
graduated in centimeters

1/2



# Forceps For Obstetrics



**04-6385**  
HERFF  
membrane puncturing  
forceps  
255 mm



**04-6386**  
IOWA  
membrane puncturing  
forceps  
260 mm




# Scissors of Obstetrics, Umbilical cord clamp



**04-6387**  
SCHUMACHER  
Umbilical cord  
scissors  
155 mm




**04-6388**  
Umbilical cord  
scissors  
105 mm




**04-6389 - 04-6390**  
BRAUN-STADLER  
Episiotomy Scissors

**04-6389**  
145 mm



**04-6390**  
200 mm



**04-6391**  
Episiotomy Scissors  
240 mm



# Obstetrics - Obstetrical Forceps



1/3

**04-6392**  
NAEGELE  
355 mm



1/3

**04-6393**  
NAEGELE  
400 mm



1/3

**04-6394**  
SIMPSON-BRAUN  
305mm



# Obstetrics - Obstetrical Forceps



1/3

**04-6395**  
SIMPSON-BRAUN  
350mm

For caesarean section



1/3

**04-6396**  
SIMPSON  
235mm

For caesarean  
section



1/3

**04-6397**  
WRIGLEY  
250 mm



# Obstetrics Forceps




**04-6398**  
**GREEN-ARMYTAGE**  
Caesarean section haemostatic  
forceps, for edges of the transverse  
incision in the lower segment  
220 mm

**04-6399**  
**WILLETT**  
Scalp Flap forceps  
185 mm

**04-6400**  
**GAUSS**  
Scalp Flap forceps  
260 mm




# Obstetrics-Placenta and Ovum Forceps






# Obstetrics-Placenta and Ovum Forceps



04-6408  
SAENGER  
285 mm



04-6409  
SAENGER  
265 mm




04-6410 - 04-6411  
LEWKOWITZ  
245 mm





# Obstetrics-Placenta and Ovum Forceps, Perforator, Cranioclast



04-6412  
KELLY  
320mm



04-6413  
SMELLIE  
Perforator  
270mm

1 jaws fenestrated



04-6414  
BRAUN  
Cranioclast  
420mm