

Section - 07



Pages 393-432

- ★ Plaster instruments
- ★ Suction instruments
- ★ Tungsten carbide Inserts
- ★ Hospital wares



Plaster Instruments

Plaster Knives



with ebony handle

07-9001
HOPKINS
200 mm



07-9002
ESMARCH
180 mm

Plaster Saw



07-9003
ENGEL
160 mm

Bandage and First-Aid scissors



07-9004 - 07-9012 145mm

- | | | | |
|---------|-------------|---------|-------------|
| 07-9004 | black | 07-9008 | yellow |
| 07-9005 | red | 07-9009 | white |
| 07-9006 | blue | 07-9010 | orange |
| 07-9007 | light green | 07-9011 | voilet |
| | | 07-9012 | bright pink |



07-9013 - 07-9021 180mm

- | | | | |
|---------|-------------|---------|-------------|
| 07-9013 | black | 07-9017 | yellow |
| 07-9014 | red | 07-9018 | white |
| 07-9015 | blue | 07-9019 | orange |
| 07-9016 | light green | 07-9020 | voilet |
| | | 07-9021 | bright pink |




Plaster Scissors

For small plaster dressings



07-9022
160 mm



07-9023 - 07-9024
ESMARCH
07-9023 200mm
07-9024 230mm



07-9025
BERGMANN
230 mm



07-9026
BRUNS
235 mm



07-9027
235 mm

also to be used as
bandage scissors



07-9028
SEUTIN
230 mm


also to be used as
bandage scissors



Plaster Shears

Special features:

The construction of the open jaws avoids tiresome clogging by plaster remnants and allows furthermore an easy cleaning. The shears do not glide off the plaster owing to the deep teething on the jaws. The knives with hollow grinding have a better cutting ability and therefore less force is required. The knives are guided by ball bearings.



07-9029
STILLE
260 mm



07-9030
STILLE
360 mm



STILLE Plaster Shears

Blade with reinforced back



1/2

07-9032
STILLE
230 mm

07-9033
Spare blade



1/2

07-9034
STILLE
370 mm

07-9035
Spare blade



1/2

07-9036
STILLE
240 mm

07-9037
Spare blade



Plaster Spreader, Plaster Breakers



1/2

07-9038
HENNIG
270 mm



1/2

07-9039
MOD.USA
230 mm



1/2

07-9040
290 mm



with removable jaws



1/2

07-9041 - 07-9042
WOLFF
07-9041 180 mm
07-9042 240 mm
plaster breakers