



# Suture Instruments - Reverdin Needles




02-3909

02-3910

02-3911

02-3912

02-3909 - 02-3912  
REVERDIN  
140 mm



02-3913

02-3914

02-3915

02-3916

02-3917

02-3918

02-3919

02-3920

02-3913 - 02-3916  
REVERDIN  
195 mm

02-3917 - 02-3920  
REVERDIN  
220 mm



# Suture Instruments - Ligature Needles







# Suture Instruments - Ligature Needles



1/2  
**02-3933**  
HURD  
for right hand  
sharp  
210 mm



1/2  
**02-3934**  
HURD  
for left hand  
sharp  
210 mm



1/2  
**02-3935**  
DESCHAMPS  
for right hand  
sharp  
215 mm



1/2  
**02-3936**  
DESCHAMPS  
for left hand  
sharp  
215 mm




1/2  
**02-3937**  
JONESCO  
Needle for wire suture  
185 mm



1/2  
**02-3938**  
CLAUS VAN EICKEN  
Needle for wire suture  
250 mm




1/2  
**02-3939**  
150 mm



1/2  
**02-3940**  
SYME  
170 mm




1/2  
**02-3941**  
COOPER  
205 mm



1/2  
**02-3942**  
KOCHER  
205 mm



# Ligature Needles



**02-3943**  
BRUNNER  
240 mm

**02-3944**  
SCHMIEDEN-DICK  
Malleable / Flexible  
320 mm

**02-3945 - 02-3946**  
BRUNNER  
Guiding Probes  
300 mm

**02-3948**

**02-3947**


**02-3945**

**02-3946**






# Guiding Probes and Suture Instruments





# Directors and Probes




**02-3962 - 02-3964**  
 02-3962 130 mm  
 02-3963 145 mm  
 02-3964 160 mm  
 Directors

**02-3965**  
 DOYEN  
 145 mm  
 Directors

**02-3966**  
 NELATON  
 160 mm  
 Directors

**02-3967**  
 Director with tip  
 Diam. Of tip 3,3 mm  
 140 mm

**02-3968**  
 BRODIE  
 Director with tip,  
 Diam. Of tip 3,5 mm  
 200 mm



**02-3969**  
 LOCKHART-MUMMERY  
 165 mm  
 Straight

**02-3970**  
 LOCKHART-MUMMERY  
 165 mm  
 45° Curved

**02-3971**  
 LOCKHART-MUMMERY  
 165 mm  
 90° Curved

**02-3972**  
 LOCKHART-MUMMERY  
 165 mm  
 Retrograde Curved

**02-3973**  
 Director with wire  
 loop, malleable  
 190 mm





# Probes

02-3974  
2,5 mm  $\varnothing$   
02-3975  
2,5 mm  $\varnothing$   
02-3976  
3,0 mm  $\varnothing$



1/2

02-3974 - 02-3976  
02-3974 145 mm  
02-3975 160 mm  
02-3976 200 mm  
Probes with eye

1,5 mm  $\varnothing$



1/2

02-3977  
Myrtle leaf probe  
145 mm

2,0 mm  $\varnothing$



1/2

02-3978 - 02-3982  
Myrtle leaf probe  
02-3978 130 mm  
02-3979 145 mm  
02-3980 160 mm  
02-3981 180 mm  
02-3982 200 mm

diam. of tips			Length
1 mm	1,5 mm	2 mm	
-	02-3984	-	160mm
02-3985	02-3986	02-3991	145mm
-	02-3987	02-3992	160mm
-	02-3988	02-3993	160mm
-	02-3989	02-3994	160mm
-	02-3990	-	160mm
-	-	02-3995	160mm

02-3983  
Probes



1/2

02-3996  
Double ended probe,  
malleable, diam. of tips  
3 mm with hole  
190 mm



02-3997  
250mm 300mm  
Double ended probes,  
tin, malleable, diam. of tips  
3,5 mm with hole



1/4

02-3998  
Double ended probe,  
tin, malleable,  
diam. of tips 4 mm  
400 mm



1/3

$\varnothing$  tip  
02-3999 - 02-4001  
02-3999 1,5 mm  
02-4000 2,0 mm  
02-4001 2,5 mm  
300 mm Probes

















# Probes - Spatulas - Cotton Applicators

<p>1/2</p>	<p>1/2</p>	<p>1/2</p>	<p>1/2</p>	<p>1/2</p> <p>German Silver</p>	
<p><b>02-4002</b>                  STACKE                  Silver                  100 mm                  ∅ 0.8mm</p>	<p><b>02-4003</b>                  KILLIAN                  165 mm                  ∅ 1.5mm</p>	<p><b>02-4004</b>                  160 mm                  ∅ 1.4mm</p>	<p><b>02-4005</b>                  SEIFFERT                  Silver                  180 mm                  ∅ 0.5mm</p>	<p><b>02-4006</b>                  BOWMAN                  130 mm</p>	
<p><b>02-4007 - 02-4011</b>                  02-4007 0.7+0.8mm ∅ Fig. 00/0                  02-4008 0.9+1.1mm ∅ Fig. 1/2                  02-4009 1.3+1.4mm ∅ Fig. 3/4                  02-4010 1.5+1.6mm ∅ Fig. 5/6                  02-4011 1.8+1.9mm ∅ Fig. 7/8</p>					
<p><b>02-4012</b>                  SPATULA                  145 mm</p>	<p><b>02-4013</b>                  SPATULA                  145 mm</p>	<p><b>02-4014</b>                  Modell PICCOLO                  Round Serrated end                  110 mm 0.7mm</p>	<p><b>02-4015</b>                  Modell PICCOLO                  Triangular Serrated end                  110 mm 0.7mm</p>	<p><b>02-4016</b>                  FARRELL                  Round Serrated end                  110 mm 1.2mm</p>	<p><b>02-4017</b>                  LATHBURY-FARRELL                  Triangular Serrated end                  110 mm 1.4mm</p>





# Cotton Applicators

 <p><b>02-4018</b> BUCK Rough, spiral end ∅120 mm 0.9mm</p>	 <p><b>02-4019</b> Triangular serrated end ∅140 mm 0.9mm</p>	 <p><b>02-4020</b> Rough, spiral end ∅140 mm 1.0mm</p>	 <p><b>02-4021</b> Rough, spiral end ∅140 mm 1.2mm</p>	 <p><b>02-4022</b> Rough, spiral end ∅140 mm 1.8mm</p>	 <p><b>02-4023</b> FARRELL Round, serrated end ∅140 mm 1.2mm</p>	 <p><b>02-4024</b> LATHBURY-FARRELL Triangular, serrated end 140 mm ∅1.4mm</p> 
 <p><b>02-4025</b> FARRELL Rough, spiral end ∅160 mm 1.2mm</p>	 <p><b>02-4026</b> LATHBURY-FARRELL Triangular, serrated end 160 mm ∅1.2mm</p>	 <p><b>02-4027</b> Triangular, serrated end ∅160 mm 1.4mm</p>	 <p><b>02-4028</b> Rough, spiral end ∅180 mm 1.5mm</p>	 <p><b>02-4029</b> Rough, spiral end ∅180 mm 1.8mm</p>	 <p><b>02-4030</b> Rough, spiral end ∅270 mm 2.5mm</p>	



# Clip Applying Forceps, Clip Applying - and Removing Forceps



½

**02-4031**  
MICHEL  
125 mm



½

**02-4032**  
HEGENBARTH  
120 mm



½

**02-4033**  
WACHENFELDT  
120 mm



½

**02-4034**  
WACHENFELDT  
120 mm



½

**02-4035**  
MICHEL  
125 mm



½

**02-4036**  
MICHEL modif.  
130 mm



½

**02-4037**  
HEATH  
135 mm



½

**02-4038**  
COLLIN  
135 mm





## Suturing Instrument Set, Pocket Instrument Set



### 02-4039

Suturing instruments set, in metal case consisting of:  
 1 michel clip applying and removing forceps,  
 1 Scissors, curved, sharp/sharp, 115mm  
 25 MICHEL clips, 16 x 3mm 1Delicate tissue forceps, 130mm  
 1 Scalpel, bellied 1 MATHIEU needle holder, 140mm  
 1 Perforated metal case, stainless  
 containing:  
 1 Suture needle  
 1 Suture needle  
 1 Suture needle



### 02-4040


Pocket instrument set, in metal case consisting of:

#### Bottom tray

1 MATHIEU needle holder, 140mm  
 1 MICHEL clip forceps, combined  
 1 Myrtle leaf probe, 130mm  
 1 Sharp double curette  
 1 Skin hook  
 1 Detachable scalpel handle, no.4  
 1 Perforated metal case, containing  
 6 Scalpel blades  
 1 Perforated metal case, stainless, containing  
 1 Suture needle  
 1 Suture needle  
 1 Suture needle

#### Top tray

1 Tissue forceps, 130mm  
 1 KOCHER haemostatic forceps, 130mm  
 1 Straight surgical scissors, 130mm  
 1 Dissecting forceps, 130mm  
 1 FEILCHENFELD splinter forceps, 115mm  
 25 MICHEL clip, 16 x 3mm  
 1 Flacon of suture skin, sterile in alcohol



### 02-4041

Pocket instrument set, in metal case consisting of:

#### Bottom tray

1 MAYO-HEGER needle holder, 150mm  
 1 MICHEL clip forceps, combined 1 Myrtle leaf probe, 130mm  
 1 Sharp double curette 1 Skin hook 1 Detachable scalpel handle, no. 4  
 1 Perforated metal case, consisting  
 6 Scalpel blades  
 1 Perforated metal case, stainless, consisting:  
 1 Suture needle  
 1 Suture needle  
 1 Suture needle

#### Top tray


1 Tissue forceps, 130mm  
 1 CRILE haemostatic forceps, 140mm  
 1 Straight surgical scissors, 130mm  
 1 Dissecting forceps, 130mm  
 1 FEILCHENFELD splinter forceps, 115mm  
 25 MICHEL clip, 16 x 3mm  
 1 Flacon of suture silk, sterile in alcohol







## Pocket Instrument Set, Venae Sectio Set



**02-4042**

Pocket instrument set in metal case 170 x 90mm  
canvas roll only, sterilizable in autoclave.  
consisting of:

- 1 Myrtle leaf probe, 130mm
- 1 Sharp double curette
- 1 FEILCHENFELD splinter forceps, 115mm
- 1 Dissecting forceps, 130mm
- 1 Tissue forceps, 130mm

- 1 Straight surgical scissors, 130mm
- 1 Scalpel handle no. 4
- 1 Perforated metal case, consisting  
6 scalpel blades
- 1 KOCHER haemostatic forceps, 130mm
- 1 Skin hook
- 1 MATHIEU needle holder, 140mm
- 1 Tube of suture skin, sterile in alcohol
- 1 MICHEL clip forceps, 120mm




- 25 MICHEL clip, 16 x 3mm
- 1 Perforated metal case, stainless,  
consisting:
- 1 Suture needle
- 1 Suture needle
- 1 Suture needle

**02-4043**

Canvas roll only, sterilizable in  
autoclave


Suitable Metal Cases



**02-4044**

Venae sectio set, canvas roll only, sterilizable in autoclave.  
consisting of:

- 2 Tissue forceps, 1 x 2 teeth, 140mm
- 1 Dissecting forceps, 145mm
- 1 FEILCHENFELD splinter forceps, 115mm
- 1 Straight scissors, sharp, 110mm
- 1 Curved scissors, blunt, 115mm
- 3 HALSTED mosquito forceps, curved on flat, 125mm



250 x 125 x 55mm

233 x 148 x 57mm

- 1 DESCHAMPS needle, blunt, 210mm
- 1 DOYEN Director, curved, 140mm
- 2 Wound retractors, 4 blunt prongs, 160mm
- 2 DESMARRES wound retractors, 135mm
- 1 FINSSEN retractors, 2 x3 blunt prongs, 80mm
- 1 Scalpel handle no. 4
- 4 Sterile scalpel blades, no 18, single peel pack
- 1 Needle box, perforated, round, diam. 65mm
- 1 HEGAR-MAYO DUROGRIP needle holder, 180mm


**02-4045**

Canvas roll only, sterilizable in autoclave





## Venae Sectio Set



### 02-4046

Venae sectio set, in metal case  
consisting of:

#### Bottom tray

- 1 MATHIEU needle holder, fine pattern 140mm
- 1 Ligature needle, fine short pattern
- 1 DOYEN director, curved
- 1 HEGENBARTH clip forceps,
- 1 ELSCHNIG tissue forceps, 1 x 2 teeth fine, 105mm
- 1 Dressing forceps, fine 105mm
- 1 Tissue forceps, 1 x 2 teeth, 110mm
- 1 Vessel scissors, straight, with 1 ball point, 110mm
- 1 Vessel scissors, curved, pointed 105mm
- 1 FINSSEN retractor
- 1 Perforated needle case

#### Top tray

- 1 Scalpel, pointed
- 1 Scalpel, bellied
- 2 HALSTED mosquito forceps, curved, 125mm
- 2 HALSTED mosquito forceps, straight, 1x 2 teeth, 125mm
- 3 BACKHAUS towel forceps, 90mm
- 2 Flacon holder